


Titanium® Ti-6Al-4V Alloy Zimmer MotionLoc Screws			Collar Section			Shaft Section				Front Section		
Screw (Sterile and Non Sterile)	Description	Pilot Hole Size	Collar Diameter	Collar Length	Collar Diameter/ Pilot Hole Clearance per side	Trailing Thread Major Diameter (central thread)	Trailing Thread Minor Diameter (central thread)	Trailing Thread Length	Diameter of Undercut Between Collar and Trailing Thread	Leading Thread Major Diameter	Leading Thread Minor Diameter	Leading Thread Length (front section)
02.02162.024 and 02.03162.024	4.0mm X 24mm Ti Zimmer MotionLoc Screw	3.3mm max size (LMC) 3.282mm min size (MMC)	2.6mm max (MMC) 2.45mm min (LMC)	6mm	.425mm max .341mm min	4mm max 3.85mm min	2.6mm max 2.45mm min	5mm	N/A	4mm max 3.85mm min	3.4mm max 3.25mm min	8mm
02.02162.026 and 02.03162.026	4.0mm X 26mm Ti Zimmer MotionLoc Screw			7mm				6mm				8mm
02.02162.028 and 02.03162.028	4.0mm X 28mm Ti Zimmer MotionLoc Screw			7mm				7mm				9mm
02.02162.030 and 02.03162.030	4.0mm X 30mm Ti Zimmer MotionLoc Screw			7mm				8mm				10mm
02.02162.032 and 02.03162.032	4.0mm X 32mm Ti Zimmer MotionLoc Screw			7mm				10mm				10mm
02.02162.034 and 02.03162.034	4.0mm X 34mm Ti Zimmer MotionLoc Screw			8mm				11mm				10mm
02.02162.036 and 02.03162.036	4.0mm X 36mm Ti Zimmer MotionLoc Screw			9mm				12mm				10mm
02.02162.038 and 02.03162.038	4.0mm X 38mm Ti Zimmer MotionLoc Screw			9mm				11mm				12mm
02.02162.040 and 02.03162.040	4.0mm X 40mm Ti Zimmer MotionLoc Screw			9mm				13mm				12mm
02.02162.042 and 02.03162.042	4.0mm X 42mm Ti Zimmer MotionLoc Screw			9mm				15mm				12mm
02.02162.044 and 02.03162.044	4.0mm X 44mm Ti Zimmer MotionLoc Screw			9mm				17mm				12mm
02.02162.046 and 02.03162.046	4.0mm X 46mm Ti Zimmer MotionLoc Screw			9mm				19mm				12mm



+H124973161003001/\$101018J10Y

97-3161-003-00 1008-T15 1ML Printed in USA ©2010 Zimmer, Inc.



Titanium <i>Zimmer MotionLoc</i> Screws			Collar Section			Shaft Section				Front Section		
Screw (Sterile and Non Sterile)	Description	Pilot Hole Size	Collar Diameter	Collar Length	Collar Diameter/ Pilot Hole Clearance per side	Trailing Thread Major Diameter (central thread)	Trailing Thread Minor Diameter (central thread)	Trailing Thread Length	Diameter of Undercut Between Collar and Trailing Thread	Leading Thread Major Diameter	Leading Thread Minor Diameter	Leading Thread Length (front section)
02.02161.030 and 02.03161.030	5.0mm X 30mm Ti <i>Zimmer</i> <i>MotionLoc</i> Screw	4.3mm max size (LMC) 4.282mm min (MMC)	3.45mm max (MMC) 3.35mm min (LMC)	5mm	.475mm max .416mm min	5mm max 4.85mm min	3.2mm max 3.05mm min	11mm	3.275mm max 3.225mm min	5mm max 4.85mm min	4.4mm max 4.25mm min	10mm
02.02161.032 and 02.03161.032	5.0mm X 32mm Ti <i>Zimmer</i> <i>MotionLoc</i> Screw			6mm				12mm				10mm
02.02161.034 and 02.03161.034	5.0mm X 34mm Ti <i>Zimmer</i> <i>MotionLoc</i> Screw			7mm				13mm				10mm
02.02161.036 and 02.03161.036	5.0mm X 36mm Ti <i>Zimmer</i> <i>MotionLoc</i> Screw			8mm				14mm				10mm
02.02161.038 and 02.03161.038	5.0mm X 38mm Ti <i>Zimmer</i> <i>MotionLoc</i> Screw			8mm				13mm				13mm
02.02161.040 and 02.03161.040	5.0mm X 40mm Ti <i>Zimmer</i> <i>MotionLoc</i> Screw			8mm				15mm				13mm
02.02161.042 and 02.03161.042	5.0mm X 42mm Ti <i>Zimmer</i> <i>MotionLoc</i> Screw			8mm				17mm				13mm
02.02161.044 and 02.03161.044	5.0mm X 44mm Ti <i>Zimmer</i> <i>MotionLoc</i> Screw			8mm				19mm				13mm
02.02161.046 and 02.03161.046	5.0mm X 46mm Ti <i>Zimmer</i> <i>MotionLoc</i> Screw			8mm				21mm				13mm
02.02161.048 and 02.03161.048	5.0mm X 48mm Ti <i>Zimmer</i> <i>MotionLoc</i> Screw			8mm				23mm				13mm
02.02161.050 and 02.03161.050	5.0mm X 50mm Ti <i>Zimmer</i> <i>MotionLoc</i> Screw			8mm				25mm				13mm
02.02161.052 and 02.03161.052	5.0mm X 52mm Ti <i>Zimmer</i> <i>MotionLoc</i> Screw			8mm				27mm				13mm
02.02161.054 and 02.03161.054	5.0mm X 54mm Ti <i>Zimmer</i> <i>MotionLoc</i> Screw			8mm				29mm				13mm
02.02161.056 and 02.03161.056	5.0mm X 56mm Ti <i>Zimmer</i> <i>MotionLoc</i> Screw			8mm				31mm				13mm
02.02161.058 and 02.03161.058	5.0mm X 58mm Ti <i>Zimmer</i> <i>MotionLoc</i> Screw			8mm				33mm				13mm
02.02161.060 and 02.03161.060	5.0mm X 60mm Ti <i>Zimmer</i> <i>MotionLoc</i> Screw			8mm				35mm				13mm